

PORTABLE AUTOCLAVE



YXQ-LS-18SII



YXQ-LS-24SII

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

PORTABLE AUTOCLAVE >>

YXQ-LS-18SII/24SII

Product Features >

- Microcomputer automatic control;
- Highlight Nixie tube display ;
- Electric heating mode ;
- Automatic control of sterilization process ;
- SUS304 material ;
- Water level safety float ;
- Self expanding sealing ring are used for the inner tank.



YXQ-LS-18SII



YXQ-LS-24SII

Model	YXQ-LS-18SII	YXQ-LS-24SII
Volume [L]	18	24
Power [kW]	2	2
Power Supply [V, Hz]	220±10%, 50±2%	220±10%, 50±2%
Max. Working / Design Temp. [°C]	126 / 128	126 / 128
Max. Working / Design Pressure [MPa]	0.142 / 0.165	0.142 / 0.165
Temperature Resolution[°C]	0.1	0.1
Timer Range [min]	0~999	0~999
Internal Dimensions [mm]	Ø300×330	Ø300×430
Housing Dimensions [mm]	Ø325×555	Ø325×660
Basket Dimensions [mm]	–	–
Number of Baskets [pc]	–	–
Net Weight [kg]	14	16
Gross Weight [kg]	16	18
Package Size [mm]	420×450×560	420×450×680