



PORTABLE AUTOCLAVE



Shanghai Boxun Medical Biological Instrument Corp

Add:No.599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel:021-56980111 Fax:021-56303876

E-mail:export@shbxyl.com

web:www.sh-bx.com



www.sh-bx.com

PORTABLE AUTOCLAVE >>

YXQ-LS-18SI

Product Features >

- Electric heating type;
- Self-control when lack of water;
- Self-expanding type seal ring;
- Self-exhaust when over pressure;
- Temperature adjustable (110°C-126°C);
- Sterilizing time adjustable (0-60min).



Pressure Gage



Water Level Safety Float



YXQ-LS-18SI

Model	YXQ-LS-18SI
Volume [L]	18
Power [kW]	2
Power Supply [V, Hz]	220±10%, 50±2%
Max. Working / Design Temp. [°C]	126 / 128
Max. Working / Design Pressure [MPa]	0.142 / 0.165
Timer Range [min]	–
Temperature Resolution[°C]	0~60
Internal Dimensions [mm]	Ø280×290
Housing Dimensions [mm]	Ø390×530
Basket Dimensions [mm]	–
Number of Baskets [pc]	–
Net Weight [kg]	19
Gross Weight [kg]	21
Package Size [mm]	450×450×640