

VERTICAL AUTOCLAVE



Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

VERTICAL AUTOCLAVE >>

BXM-30R

Product Features >

- Made of SUS304;
- Microcomputer intelligent control;
- Self-expanding type seal ring;
- Exhaust cold air automatically;
- Self-control when lack of water;
- Over temperature protection;
- Self-exhaust when over pressure;
- Buzzer reminds when sterilization ends.



Basket



Stainless Steel Bucket

Model	BXM-30R
Volume [L]	29
Power [kW]	2
Power Supply [V, Hz]	220±10%, 50±2%
Max. Working / Design Temp. [°C]	126 / 128
Max. Working / Design Pressure [MPa]	0.142 / 0.165
Temperature Resolution[°C]	1
Timer Range [min]	0~120
Internal Dimensions [mm]	Ø300×420
Housing Dimensions [mm]	395×395×780
Basket Dimensions [mm]	Ø280×340
Number of Baskets [pc]	1
Net Weight [kg]	35
Gross Weight [kg]	43
Package Size [mm]	540×540×950