



MICROBIOLOGICAL INCUBATOR



BXP-S Series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com




www.sh-bx.com

MICROBIOLOGICAL INCUBATOR >>

BXP-S Series

Product Features

- 5 inch touch screen, BRIGHT II control system, automatically compensate parameters based on the environment;
- Temperature record data and curves, alarm recording information can be viewed in real time;
- Multi-level password management function;
- With reservation function;
- Store data for more than 10 years, and export them by U disk;
- Double overtemperature alarm (a. controlled by temperature sensor system, b. independent temperature control switch);
- The 6-stage speed regulation of the fan;
- 10 adjustable exhaust rates;
- Preheating chamber design, to ensure rapid heating and good heat distribution effect;
- 5cm fully enclosed aluminum foil insulation layer, suitable for use in clean laboratories;
- Double door design for use in clean laboratories;
- Equipped with 2 RS-485 interfaces and 1 USB interface;
- Optional: Printer, online monitoring platform, inlet and outlet air filter.



Reservation Timing Interface



Operating Main Interface



Data Curve Viewing



USB Export Interface

Model	BXP-65S	BXP-130S	BXP-280S	BXP-450S	BXP-530S
Volume [L]	65	130	280	450	530
Temp. Range [°C]	RT+5~100				
Temp. Fluctuation at 37°C [°C]	±0.2				
Temp. Uniformity at 37°C [°C]	±0.6				
Heating-up Time to 60°C [min]	16	20	15	25	25
Recovery time after 30 seconds door open at 60°C [min]	14	13	10	11	11
Power [W]	550	550	1000	1500	1000
Internal Dimensions W×D×H [mm]	410×350×450	570×400×570	720×530×745	600×580×1300	1070×595×835
Housing Dimensions W×D×H [mm]	575×615×635	735×665×755	885×795×930	774×845×1574	1235×850×1080
Number of Shelves (std./max.) [pc]	2/5	2/7	3/9	4/13	3/12
Net Weight/Gross Weight [kg]	52/73	67/88	102/125	150/179	167/190
Package Size [mm]	695*768*787	855*818*907	1005*948*1082	894*1005*1724	854*965*1614