

VACUUM OVEN



DZF Series

Shanghai Boxun Medical Biological Instrument Corp
上海博迅医疗生物仪器股份有限公司

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

Vacuum Oven

DZF Series

Product Feature

- The shell is made of cold rolled steel plate, the surface is electrostatically sprayed, the inner is made of stainless steel;
- Double screen high brightness Nixie tube display, touch buttons ;
- Microcomputer single-chip temperature control system, temperature control, timing, over temperature alarm functions;
- Timing and timing functions ;
- The door can be adjusted, a silicone rubber door seal, toughened, bulletproof double glass door;
- Power supply recovery function, data will not be lost when power failure and system halted;
- Vacuum pump is optional.

Model	DZF-6020MBE	DZF-6050MBE
Volume [L]	25	53
Temp. Range [°C]	RT+5~250	RT+5~250
Temp. Fluctuation at 105°C [°C]	±1	±1
Vacuum Degree [Pa]	133	133
Power [W]	900	1400
Power Supply [V, Hz]	220±10, 50/60	220±10, 50/60
Timer Range [min]	0~999	0~999
Internal Dimensions W×D×H [mm]	300×296×275	415×370×344
Housing Dimensions W×D×H [mm]	590×430×475	710×505×550
Number of Shelves [pc]	2	2
Net Weight [kg]	40	64
Gross Weight [kg]	49	74
Package Size [mm]	690*570*630	810*650*700